DynaCox® evolution

Hüftorthese

Hip Orthosis

Patienteninformation

Patient Information

CE

www.wjt-ortho.com

WJT
ORTHO

◆ Liebe Anwenderin, lieber Anwender,

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◆ Dear Patient,

You have received an artificial hip joint. This kind of surgery is a safe standard procedure which is performed more than 280.000 times* a year in Germany.

By implanting an artificial hip joint, pain that was previously triggered by motion can in many cases be completely eliminated or at least significantly reduced. Your new hip joint is not the end in a chain of therapeutic actions; it is rather the beginning of new opportunities for activity!

Aside from giving instructions on how to use the DynaCox® evolution, this brochure is intended to answer some questions about your artificial hip joint and give you some useful advice for everyday life with it. The brochure contains extracts of the publication "Die Endoprothesen-Schule: Ratgeber für Patienten" (Endoprosthesis Training: A Patient's Guide) by Prof. Dr. med. J. Jerosch and Prof. Dr. med. J. Heisel, which we recommend for further reading (ISBN: 3-934849-84-9). It is available through your local bookshop.

Below, extracts from the publication are indicated by italic type.

◆ Contents

Description	EN-3
• Basic Rules	EN-3
• Artificial Joints and Stress	EN-4
• After Surgery	EN-4
• The 10 Most Important Rules of Conduct for Life with an Artificial Joint	EN-5
• Everyday Activities	EN-6
• Sleeping	EN-6
• Sitting	EN-6
• Standing	EN-7
• Walking	EN-7
• Walking on Stairs with Crutches	EN-7
• Lifting Objects	EN-8
• Carrying Loads	EN-8
• Getting Dressed and Undressed	EN-8
• Shoes	EN-8
Indication	EN-8
Contraindication	EN-9
Function	EN-9
Medical Device	EN-9
Safety	EN-9
Application/Adjustment	EN-10
Comfort/Hygiene	EN-11
Material	EN-11
Product Care/Cleaning/Maintenance	EN-12
Storage/Disposal	EN-12
Legend to Symbols	EN-12

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◆ Description

The DynaCox® evolution hip brace is designed to protect your hip and minimize the risk of a dislocation.

The orthosis limits how much you can move, so movements which would stimulate a luxation are avoided. Furthermore, the dynamic functionality of the orthosis helps you maintain a natural gait and promotes an active use of the musculature.

The brace acts as a reminder when you try to move in a way that might cause a luxation.

Do not test the brace to see how far you can bend your hip while wearing the orthosis. Also, do not try to cross your legs against the pressure applied by the orthosis.

In addition to following your attending physician's advice, it will be necessary to follow some general rules. This brochure will provide you with information on how to avoid luxation and how to safely use the DynaCox® evolution. Please pay heed to „The 10 Most Important Rules for Life with an Artificial Joint“.

• Basic Rules

- Before using this device, an individual adaptation of the DynaCox® evolution and training in its use by an orthopedic technician according to the medical instructions is essential.
- If an inappropriate size is selected, if the individual adjustment of the brace is neglected or if monitoring is insufficient, the success of the therapy may be strongly compromised (possibly resulting in a total lack of luxation prophylaxis).
- Check the orthosis for damage (e.g. loose connections, fasteners, belts, etc) each time it is put on.
- Go to see your doctor or orthopedic technician immediately if you notice performance deficiencies (e. g. feeling that the orthosis does not stabilize the hip as much as usual) or the DynaCox® evolution has been subjected to a disproportionately high stress.
- Have the orthosis checked for correct fit and function by an orthopedist or orthopaedic technician at regular intervals. The intervals between checks have to be determined by an orthopedist or orthopedic technician.

• Application Periods of the DynaCox® evolution Orthosis

Your attending physician will talk to you about the estimated application period and the necessary daily wearing time for the DynaCox® evolution. For optimal luxation protection it is imperative to observe your doctor's instructions. Depending on how high the risk of a luxation is, it may be necessary to wear the DynaCox® evolution around the clock for a defined period of time.

• Artificial Joints and Stress

What stresses is the artificial joint subjected to?

Due to its carrying function as a locomotion organ the hip joint is subjected to highly variable (...) stresses (...) induced by the weight of the body. (...)

To evaluate the forces acting on the joint during everyday activities (...) the stresses that arise during walking and running movements are crucial. Over the course of a human life, these loads add up (...) to appr.1-2.5 million load changes, i.e. the load shifting from one leg to the other.

During stance phase, which is the phase in walking during which most of the weight rests on one leg, the hip is loaded with 2.5 to 3.5 times the body weight. During activities such as jogging, this weight can increase to up to six times the body weight.

Why does an artificial joint wear out?

Despite all the progress in technology and surgical medicine, it may not be forgotten that an artificial joint is a technical component which is subject to wear in the course of its service life.

Usually, a period of 10 to 15 years is given in this context. (...) One of the main problems is the direct contact of joint surfaces, which are subject to constant friction.

Special tribological pairings are used to keep the phenomenon of friction to a minimum.

However, in the course of constant progress in implant technology, a further increase of the average performance life can be expected in the future.

• After Surgery

What needs to be taken into consideration?

Soon after surgery, the time after your stay in hospital should be planned in detail. Discuss your domestic environment with a physiotherapist or ergotherapist while still in hospital. Also consider the various kinds of secondary treatment available (rehabilitation measures). Insist on receiving the secondary treatment concepts in writing before leaving the hospital. Do not leave the hospital before you have had a detailed concluding

meeting with your operating physician and your physiotherapist.

Do not forget to always carry your endoprosthesis-passport on you. In case of emergency (accident, luxations) it will provide crucial information to the attending doctor.

Once you have been discharged from the hospital, wearing your anti-thrombosis stockings until you can fully put weight on the leg again.

In case you suddenly encounter pain in the calf or swelling of the lower leg or foot, please see your doctor immediately (danger of thrombosis!).

• The 10 Most Important Rules of Conduct for Life with an Artificial Joint

- 1** An endoprosthesis can never fully replace a natural joint!
- 2** A few weeks after surgery all normal coordinated movements are possible – but extreme movements must be avoided!
- 3** When sitting, your knee joints must not be in a position higher than your hips (danger of luxation!).
- 4** In everyday life you should load the operated leg as evenly as possible, avoid movements requiring intense strength!
- 5** Avoid carrying loads weighing more than 20% of your own body weight!
- 6** Be especially careful when there is an increased risk of falling (wet or icy conditions)!
- 7** The artificial joint must always be protected from suppurative inflammation which constitutes a dangerous complication. Therefore, an early broadspectrum antibiotic prophylaxis is needed in case dental treatments or medical measures involving the urinary tract are required!
- 8** In case of vague or, particularly, increasing pain in the area of the artificial joint, especially when during physical exertion you should immediately see your attending physician!
- 9** Even if you do not observe any significant discomfort or ailments, you should have your artificial joint examined once a year by your physician!
- 10** Carry your endoprosthesis-passport and your emergency health card with you at all times (especially on flights due to the safety checks at airports)!

• Everyday Activities

In your daily life you move your body in many different ways.

Some are helpful for building your muscles and staying mobile; others should be avoided with regard to maintaining the service life of the artificial joint for as long as possible and to preventing complications such as luxations.

The advice below should be categorically followed to for the first months at home

- No hip flexion beyond 90°
- No abduction of the operated leg beyond the body's middle line
- No pronounced rotation of the hip joint to the inside or outside

In case of a timely recovery without complications (...) you may gradually increase the level of your everyday activities.

About five weeks after surgery, you may return to the often more comfortable lateral sleeping position however, when lying on the non-operated leg you should still use the abduction pillow temporarily.

Usually, lying on the operated side will be possible again about seven weeks after surgery without any problems.

As soon as you are able to flex your hip 90°, you are allowed to sit in a normal chair (starting around the 7th week). However, particularly low or soft

upholstered furniture should not yet be used. Always use your non-operated leg when getting up from a low sitting position.

*You may now also begin to put on shoes and socks yourself. Helpful items for this purpose are **long shoehorns, elastic shoelaces and dressing aids for putting on socks.***

• Sleeping

At first you should always sleep on your back with your legs slightly abducted. An abduction pillow will help you to maintain a slightly abducted leg position.

Avoid turning over in bed and all movements involving turning your body or hip.

• Sitting

After having your hip joint replaced without complications you may usually sit up beginning on the 5th day after surgery. You should refrain from flexing the hip too much (as is required for sitting on chairs of normal height).

*For this purpose an elevated chair or a sloped **arthrodesis cushion** which reduces the hip flexion angle may be used.*

Pay particular attention to the operated leg standing firmly on the floor while

sitting. If you let the leg swing freely (...), gravity will pull it over the edge of the chair even more. This will lead to an unintentional hyperextension of the leg which may cause the artificial joint to be pulled forward and out of its socket.

Avoid crossing your legs, as well as extreme inside or outside turns of the hip joint. Another useful item is an **elevated toilet seat**.

• Standing

While standing, it is important that you distribute your weight on both legs as evenly as possible. Don't turn suddenly and abruptly, but make sure to slowly turn your pelvis and legs at the same time.

• Walking

Never walk so far that you become extremely tired, but rather put in small breaks at regular intervals. Avoid leaving the house in particularly bad weather (snow, ice); if unavoidable we recommend that you use **shoe snow chains**.

• Walking on Stairs with Crutches

Learning how to use crutches for walking on stairs is one of the most important measures for you to manage everyday challenges. In the early phase after your surgery, walking upstairs will be better for you than walking downstairs.

Walking upstairs:

If there is a handrail, hold on to it with one hand and hold the crutch with your other hand. While holding on tightly with both hands, put the non-operated leg onto the first step. Pull the operated leg up onto the first step and then bring the crutch up as well.

If there is no handrail, put the foot of the non-operated leg one step higher first.

Then push off using your crutches and extend the hip and knee joint of the healthy leg. In doing so you will bring the operated leg up to the other leg on the next step, then let the crutches follow afterwards. Repeat this process for each of the following steps.

Walking downstairs:

If there is a handrail, hold on to the handrail with one hand and hold the crutch in the other hand. Now, put the crutch on the next lower step, then the operated leg and lastly the non-operated leg. If there is no handrail, first put both crutches one step lower and shift your weight to

the healthy leg at the same time. Then put the operated leg on the next lower step. Now shift your body weight to the crutches and move one step down while bending the healthy leg in the hip and the knee joint.

This step downwards is concluded by placing the healthy leg on the step where the crutches and the operated leg already are.

• Lifting Objects

Keep the operated hip joint as straight as possible. Do so by supporting yourself with one hand, e.g. on a table, and bending the knee of the non-operated leg while standing. During this process move the operated leg backwards with your hip stretched. This way your spare hand can easily reach the floor. We recommend you use a reacher, with which you will be able to pick up light objects from the floor without bending down.

• Carrying Loads

Avoid lifting and carrying loads that weigh more than 20% of your own body weight!

• Getting Dressed and Undressed

*Use **donning aids for socks and trousers** in order to avoid excessive bending.*

• Shoes

*Wear sturdy shoes without shoelaces or lace-up shoes with **elastic shoelaces** as often as possible. Use a **long shoe-horn** to put your shoes on.*

◆ Indication

- After repositioning of luxated hip joints
- Instability after T.H.R.
- Muscular imbalance
- After revision surgery

◆ Contraindication

- Heavy skin alterations

◆ Function

- Securing surgical results
- Luxation prophylaxis
- Early functional motion therapy

◆ Medical Device

The DynaCox® evolution is a medical device.

◆ ⚠ Safety

- Before this device is used, the orthosis has to be individually adapted to the patient by an experienced orthopedic technician and the patient has to be trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, application, and/or failure to monitor the correct fit of this product may cause health damage.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Do not remove the textile label from textile products.
- In case you suspect any performance malfunctions (e.g. an apparent reduction or lack of effect), immediately call on your doctor or orthopedic technician.
- All serious incidents that occur in connection with this product have to be reported

to the Wilhelm Julius Teufel GmbH and the competent authority of the member state in which you are established. A „serious incident“ is any incident that caused, may have caused, or may cause, directly or indirectly, one of the following consequences:

- a) death of a patient, user, or other persons
 - b) temporary or permanent health deterioration of a patient, user or other persons
 - c) a serious public health threat.
- Never wear the orthosis on broken skin.
 - The orthosis must be checked for proper function, correct fit and potential damage by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined individually for each patient.
 - Should the product have been subject to disproportionate stress, it must be checked for potential damage by an orthopedic technician before further use.
 - Harmful pressure caused by the orthosis may cause skin injuries or disrupt blood flow and must be avoided.
 -  This product is intended for use on one person. It may not be reused on another person.

◆ Application/Adjustment

The brace will work best if applied correctly. The brace has been adjusted in the hospital especially for you.

1. Make sure that the bandage is positioned low on your pelvis. The rigid band should be positioned well below your waist.
2. Fasten the bandage firmly around your lower abdomen, using the hand straps to help with grip. (Fig. 1)

1

3. Fasten the buckle of the rigid band. The buckle should sit well below the waist but not so low that you cannot bend your hip. Tighten the band as much as possible. (Fig. 2)

2

4. Put the clasp around your thigh and close it on the inside and outside of your leg, using the hook-and-loop strap on the back. (Fig. 3)

3

5. Stretch the bandage down your thigh to minimise creasing, and close it around your thigh. Make sure it fits as snugly and smoothly as possible. (Fig. 4)

4

6. Ensure the rigid pad is just below your hip joint and you feel it support your posture. (Fig. 5)

Your brace has now been completely applied.

The products listed below are available in your medical supply shop or through your attending orthopedic technician.

REF	Product
23.601.XXX.00	Substitution bandage for DynaCox® evolution
41.10X.X00.00	Arthrodesis cushion

Other recommended products:

Elevated toilet seat, shoe snow chains, donning aids for socks and trousers, forearm crutches, elastic shoelaces, abduction pillow, reacher.

◆ Comfort/Hygiene

Do not wear the orthosis directly on the skin, but instead over thin clothing (e.g. a t-shirt or an undershirt). Depending on the doctor's prescription it may be necessary to wear your underwear over the orthosis, in case the orthosis may not be taken off during visits to the bathroom.

◆ Material

Bandage:

48% Polyester
36% CoolMax®
16% Spandex

Orthosis component:

Steel bar
Polyamide PA 6
PE

◆ Product Care/Cleaning/ Maintenance

We recommend using a substitute bandage for washing if the daily period of use is very long.

To allow for washing, the bandage should be replaced by your specialist.

Hand wash only.

Do not iron.

Do not tumble dry.

Do not bleach.

Do not dry clean.

◆ Storage/Disposal

Store the product in its original packaging in a dry place. Protect from dust, dirt and moisture. This product does not contain any environmentally hazardous substances. It can be disposed of with regular household waste.

◆ Legend to Symbols

Manufacturer

Date of Manufacture

Distributor

Importer

EC Authorized Representative

Lot Number

Part Number

Serial Number

UDI Data Matrix

Follow User Manual

Caution

CE-Marking According to European Union Council Regulation 2017/745 Concerning Medical Devices (MDR)

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